

National Health and Nutrition Examination Survey

August 2021-August 2023 Data Documentation, Codebook, and Frequencies

Hepatitis E: IgG & IgM Antibodies (HEPE_L)

Data File: HEPE_L.xpt

First Published: February 2026

Last Revised: NA

Component Description

Hepatitis viruses constitute a major public health problem because of the morbidity and mortality associated with the acute and chronic consequences of these infections. Because of the high rate of asymptomatic infection with these viruses, information about the prevalence of these diseases is needed to monitor prevention efforts. By testing a nationally representative sample of the U.S. population, NHANES provides the most reliable estimates of age-specific prevalence needed to evaluate the effectiveness of the strategies to prevent these infections. NHANES testing for markers of infection with hepatitis viruses is used to determine secular trends in infection rates across most age and racial/ethnic groups, and provides a national picture of the epidemiologic determinants of these infections. In addition, NHANES provides the means to better define the epidemiology of other hepatitis viruses.

Hepatitis is inflammation of the liver most often caused by a virus. Viral hepatitis is a major public health problem of global importance because of the ongoing transmission of viruses that cause the disease and increased morbidity and mortality associated with the acute and chronic consequences of these infections. Global and U.S. goals have been established for elimination of viral hepatitis as a public health threat by 2030 (HHS Healthy People, 2022 and HHS 2020).

In the U.S., the most common types of viral hepatitis are hepatitis A, B, and C. Effective vaccines are available to help prevent hepatitis A and hepatitis B. No vaccine is available for hepatitis C; however, highly effective, well-tolerated treatment can cure hepatitis C virus infection. Hepatitis D virus infection is less common in the U.S. and can occur only among people with hepatitis B virus infection. Hepatitis E infection also is less common in the U.S. These five hepatitis viruses, also called hepatitides, are well-characterized for detection with laboratory assays and are monitored in U.S. public health surveillance systems.

NHANES viral hepatitis data are used to monitor progress toward goals in *Healthy People* and the Health and Human Services (HHS) *Viral Hepatitis National Strategic Plan*, which in turn support U.S. and global viral hepatitis elimination goals (HHS Healthy People, 2022 and National Academies of Sciences, Engineering, and Medicine, 2017). The viral hepatitis laboratory and interview components of NHANES complement data from outbreak, case-based surveillance, vital statistics, health care systems, and cohort studies that can provide timely, detailed, or longitudinal information for subnational geographic areas and disproportionately affected populations, such as people experiencing homelessness or living in correctional facilities; however, these sources lack information available from NHANES, such as race, ethnicity, education, income, and health status and behaviors.

Viral hepatitis data from NHANES are available beginning with NHANES II conducted in 1976-1980 for hepatitis A and hepatitis B, and with NHANES III conducted in 1988-1994 for hepatitis C, hepatitis D, and hepatitis E.

Infection with hepatitis E virus (HEV) has been responsible for large water-borne epidemics of acute disease in developing countries and for acute sporadic disease in industrialized

developed countries. In immunocompromised individuals, infection with HEV may also cause chronic infection, which may progress to fibrosis and cirrhosis.

Eligible Sample

Examined participants aged 6 years or older were eligible.

Description of Laboratory Methodology

Hepatitis E IgG Antibody (IgG Anti-HEV)

Abia HEV IgG enzyme immunoassay, which was formerly DS-EIA-ANTI-HEV-G, is an enzyme immunoassay kit intended for the detection of IgG antibodies to hepatitis E virus in human serum or plasma. During an initial incubation, HEV antibody in the sample binds with HEV antigen coated onto wells of a polystyrene stripped plate. Unbound sample is removed by washing. During a second incubation, horseradish peroxidase-labeled antibody conjugate (monoclonal mouse antibodies against human IgG) binds to any human IgG from the sample that was captured on the well. Unbound conjugate is removed by washing and a substrate solution containing tetramethylbenzidine is added to produce color. The reaction is stopped by adding a sulphuric acid solution and the optical density (OD) of each well is read. The presence or absence of IgG antibodies to hepatitis E virus is determined by the ratio of the OD of each sample to the calculated cut-off value. A sample is considered negative if the sample OD value is less than the cut-off and is considered positive if the sample OD value is greater than or equal to the cut-off.

Hepatitis E IgM Antibody (IgM Anti-HEV)

Abia HEV IgM enzyme immunoassay, which was formerly DS-EIA-ANTI-HEV-M is an enzyme immunoassay kit intended for the detection of IgM antibodies to hepatitis E virus in human serum or plasma. During an initial incubation, HEV antibody in the sample binds with HEV antigen coated onto wells of a polystyrene stripped plate. Unbound sample is removed by washing. During a second incubation, horseradish peroxidase-labeled antibody conjugate (monoclonal mouse antibodies against human IgM) binds to any human IgM from the sample that was captured on the well. Unbound conjugate is removed by washing and a substrate solution containing tetramethylbenzidine is added to produce color. The reaction is stopped by adding a sulphuric acid solution and the OD of each well is read. The presence or absence of IgM antibodies to hepatitis E virus is determined by the ratio of the OD of each sample to the calculated cut-off value. A sample is considered negative if the sample OD value is less than the cut-off and is considered positive if the sample OD value is greater than or equal to the cut-off.

Refer to the Laboratory Method Files section for a detailed description of the laboratory methods used.

There were no changes to laboratory methods, lab equipment, or lab site for this component in the August 2021-August 2023 cycle.

Laboratory Method Files

[Hepatitis E IgG Antibody](#) (February 2026)

[Hepatitis E IgM Antibody](#) (February 2026)

Laboratory Quality Assurance and Monitoring

Serum samples were processed, stored, and shipped to the Division of Viral Hepatitis, National Center for HIV/AIDS, Viral Hepatitis, STD, and TB Prevention, Centers for Disease Control and

Prevention, Atlanta, GA for analysis.

Detailed instructions on specimen collection and processing are discussed in the [NHANES Laboratory Procedures Manual \(LPM\)](#). Vials were stored under appropriate frozen (–30°C) conditions until they were shipped to Division of Viral Hepatitis, National Center for HIV/AIDS, Viral Hepatitis, STD, and TB Prevention for testing.

The NHANES quality assurance and quality control (QA/QC) protocols meet the 1988 Clinical Laboratory Improvement Act mandates. Detailed QA/QC instructions are discussed in the [NHANES LPM](#).

Mobile Examination Centers (MECs)

Laboratory team performance is monitored using several techniques. NCHS and contract consultants use a structured competency assessment evaluation during visits to evaluate both the quality of the laboratory work and the quality-control procedures. Each laboratory staff member is observed for equipment operation, specimen collection and preparation; testing procedures and constructive feedback are given to each staff member. Formal retraining sessions are conducted annually to ensure that required skill levels are maintained.

Analytical Laboratories

NHANES uses several methods to monitor the quality of the analyses performed by the contract laboratories. In the MEC, these methods include performing blind split samples collected on “dry run” sessions. In addition, contract laboratories randomly perform repeat testing on 2% of all specimens.

Data Processing and Editing

The data were reviewed. Incomplete data or improbable values were sent to the performing laboratory for confirmation.

Analytic Notes

There are over 800 laboratory tests performed on NHANES participants. However, not all participants provided biospecimens or enough volume for all the tests to be performed. The specimen availability can also vary by age or other population characteristics. Analysts should evaluate the extent of missing data in the dataset related to the outcome of interest as well as any predictor variables used in the analyses to determine whether additional re-weighting for item non-response is necessary.

Please refer to the NHANES [Analytic Guidelines](#) and the on-line NHANES [Tutorial](#) for further details on the use of sample weights and other analytic issues.

Phlebotomy Weights

For the August 2021-August 2023 cycle, analysis of nonresponse patterns for the phlebotomy component in the MEC examination revealed differences by age group and race/ethnicity, among other characteristics. For example, approximately 67% of children aged 1-17 years who were examined in the MEC provided a blood specimen through phlebotomy, while 95% of examined adults aged 18 and older provided a blood specimen. Therefore, an additional phlebotomy weight, WTPH2YR, has been included in this data release to address possible nonresponse bias. Participants who are eligible but did not provide a blood specimen have their phlebotomy weight assigned a value of “0” in their records. The phlebotomy weight should be used for analyses that use variables derived from blood analytes, and is included in all relevant data files.

Demographic and Other Related Variables

The analysis of NHANES laboratory data must be conducted using the appropriate survey design and demographic variables. The [NHANES August 2021-August 2023 Demographics File](#) contains demographic data, health indicators, and other related information collected during household interviews as well as the sample design variables. The recommended procedure for variance estimation requires use of stratum and PSU variables (SDMVSTRA and SDMVPSU, respectively) in the demographic data file.

This laboratory data file can be linked to the other NHANES data files using the unique survey participant identifier (i.e., SEQN).

Detection Limits

This data is qualitative. The use of lower limits of detection (LLODs) is not applicable.

References

- National Academies of Sciences, Engineering, and Medicine. 2017. A national strategy for the elimination of hepatitis B and C. Washington, DC: The National Academies Press. Available from: <http://www.nationalacademies.org/hmd/reports/2017/national-strategy-for-the-elimination-of-hepatitis-b-and-c.aspx>
- U.S. Department of Health and Human Services. 2020. Viral Hepatitis National Strategic Plan for the United States: A Roadmap to Elimination (2021–2025). Washington, DC. Available from: <https://www.hhs.gov/hepatitis/viral-hepatitis-national-strategic-plan/index.html>
- U.S. Department of Health and Human Services. Healthy People. 2022. Available from: <https://health.gov/our-work/national-health-initiatives/healthy-people>.

Codebook and Frequencies

SEQN - Respondent sequence number

Variable Name:	SEQN
SAS Label:	Respondent sequence number
English Text:	Respondent sequence number
Target:	Both males and females 6 YEARS - 150 YEARS

WTPH2YR - Phlebotomy 2 Year Weight

Variable Name: WTPH2YR
SAS Label: Phlebotomy 2 Year Weight
English Text: Phlebotomy 2 Year Weight
Target: Both males and females 6 YEARS - 150 YEARS

Code or Value	Value Description	Count	Cumulative	Skip to Item
4391.8220579 to 241728.85724	Range of Values	7316	7316	
0	No blood sample provided	752	8068	
.	Missing	0	8068	

LBDHEG - Hepatitis E IgG (anti-HEV)

Variable Name: LBDHEG

SAS Label: Hepatitis E IgG (anti-HEV)

English Text: Hepatitis E IgG (anti-HEV)

Target: Both males and females 6 YEARS - 150 YEARS

Code or Value	Value Description	Count	Cumulative	Skip to Item
1	Positive	836	836	
2	Negative	6284	7120	
3	Indeterminate	1	7121	
.	Missing	947	8068	

LBDHEM - Hepatitis E IgM (anti-HEV)

Variable Name: LBDHEM

SAS Label: Hepatitis E IgM (anti-HEV)

English Text: Hepatitis E IgM (anti-HEV)

Target: Both males and females 6 YEARS - 150 YEARS

Code or Value	Value Description	Count	Cumulative	Skip to Item
1	Positive	128	128	
2	Negative	708	836	
3	Indeterminate	0	836	
.	Missing	7232	8068	